

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**



## Classification

According to EN60601-1 1988, **AEROS 4600** belongs to the following classifications :

- ♦ Class I equipment
- ♦ Type B equipment
- ♦ IPX0
- ♦ Transportable equipment
- ♦ Internal power supply
- ♦ Operate continuously

## Physical specification

All specifications are approximately, maybe changed at any moment without notice.

|           |  |         |
|-----------|--|---------|
| Main unit | Height   | 280 mm  |
|           | Width  | 303 mm  |
|           | Depth  | 350 mm  |
|           | Weight   | 15 kg   |
| Cart      | Height   | 1100 mm |
|           | Width  | 650 mm  |
|           | Depth  | 650 mm  |
|           | Weight   | 23 kg   |
| Truckle   | Diameter 75 mm, the forth two truckles with arresters. |         |

## Environment requirements

|                        | Temperature | Relative humidity | Atmospheric Pressure |
|------------------------|-------------|-------------------|----------------------|
| Operation              | 5 ~ 40 °C   | ≤ 90%             | 70 ~ 106 kPa         |
| Storage/Transportation | -20 ~ 55 °C | ≤ 93%             | 50 ~ 106 kPa         |

**CAUTION:** The device should be stored at the room that is drafty and no corrosion gas exists.

**CAUTION:** When the transportation (storage) conditions are beyond the environment requirements, and the transportation (storage) state is transferred into operation state, the product only can be used after being stored in operation environment for over 24 hours.

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**

## SYSTEM TECHNICAL SPECIFICATION

### Gas supply

|             |                           |
|-------------|---------------------------|
| Type        | Compressed O <sub>2</sub> |
| Pressure    | 280 ~ 600 kPa             |
| Rating flow | max 120 l/min             |

### Power supply

|                  |  |
|------------------|--|
| Voltage          | AC 100 V ~ 240 V, 50 ~ 60 Hz DC 24V 2.3 AH |
| Input power      | ≤ 80 VA (AC 220 V, input current ≤ 0.36 A) |
| AC fuse          | 250 V 1 A Φ5X20(T) Two                     |
| DC fuse          | 250 V 6.3 A Φ5X20(T) One                   |
| Earth resistance | < 0.1 Ω                                    |

**WARNING:** If having any problem in wire protective earthing during installation and layout, the ventilator must work by internal power supply.

## VENTILATOR TECHNICAL SPECIFICATION

### Performance of ventilation

|  |  |
|--|--|
| <b>Maximum security pressure of airway system:</b> 80 cmH <sub>2</sub> O |  |
| <b>System compliance:</b>  | ≤ 4 ml/cmH <sub>2</sub> O  |
| <b>Electrical safety:</b>  | Meet requirements for Class I, Type B, equipment specified in EN60601-1 Medical Electrical equipment: Part one: General requirement for safety and EN60601-2-12 Medical electrical equipment-Part 2: Particular requirements for the safety of lung ventilators-Critical care ventilators. |
| <b>Noise of whole unit:</b>  | ≤ 65dB (A)   |
| <b>Power</b>   | Not more than 80 VA  |
| <b>Preheating time</b>   | Not more than 5 min  |
| <b>Inspiratory and expiratory resistance</b>                             | At flow of 60 l/min for adult use, inspiratory resistance ≤ 0.6 kPa; expiratory resistance ≤ 0.6 kPa<br>At flow of 30 l/min for pediatric, inspiratory resistance ≤ 0.6 kPa; expiratory resistance ≤ 0.6 kPa   |

**WARNING:** Adding accessories will add impedance when using AEROS 4600 ventilator.

### Ventilation mode

**AEROS 4600** is having the following ventilation modes:

- ◆ VCV
- ◆ VCV+Sigh
- ◆ SIMV
- ◆ SPONT
- ◆ PCV
- ◆ PSV

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**

## Setting ventilating parameters

| Item                    | Range                      | Resolution                                      | Accuracy  | Remark  |
|-------------------------|----------------------------|---|---|---|
| <b>f</b>                | 1 ~ 70 bpm                 | 1 bpm   | ±1 bpm (≤10 bpm);<br>±10% (other)   | F in VCV mode   |
| <b>V<sub>T</sub></b>    | 50 ~ 1500 mL               | 10 mL (0 ~ 1000 mL);<br>100 mL (1000 ~ 1500 mL) | ±25 mL (50 ~ 100 mL);<br>±15% (other)   | -----   |
| <b>f</b>                | 1 ~ 40 bpm                 | 1 bpm   | ±1 bpm (≤10 bpm);<br>±10% (other)   | F in SIMV mode  |
| <b>T<sub>i</sub></b>    | 0.2 ~ 6 s                  | 0.1 s   | ±0.1 s (0.2 ~ 0.5 s);<br>±20% (other)   | ----  |
| <b>T<sub>P</sub></b>    | 0 ~ 2 s                    | 0.1 s   | ±0.1 s (0 ~ 0.5 s);<br>±20% (other)   | Limited by T <sub>i</sub> , the maximum of T <sub>P</sub> is 50% T <sub>i</sub> |
| <b>FiO<sub>2</sub></b>  | 21% ~ 100%                 | 1%  | ±15%  | Invalid when no oxygen supply operated  |
| <b>P<sub>sens</sub></b> | -20 ~ 0 cmH <sub>2</sub> O | 1 cmH <sub>2</sub> O                            | ±1 cmH <sub>2</sub> O (-10 ~ 0 cmH <sub>2</sub> O);<br>±10% (other)   | -----   |
| <b>V<sub>sens</sub></b> | 2 ~ 30 l/min               | 1 l/min   | ±1 l/min (2 ~ 5 l/min);<br>±20% ( other)  | -----   |
| <b>PEEP</b>             | 0 ~ 30 cmH <sub>2</sub> O  | 1 cmH <sub>2</sub> O                            | ±2 cmH <sub>2</sub> O (≤10 cmH <sub>2</sub> O);<br>±20% (other)   | -----   |
| <b>PSV</b>              | 0 ~ 60 cmH <sub>2</sub> O  | 1 cmH <sub>2</sub> O (based on PEEP)            | The error upper limit: setting value + 6 cmH <sub>2</sub> O. The error lower limit: (setting value within 10 cmH <sub>2</sub> O) setting value – 2 cmH <sub>2</sub> O; (setting value is others) setting value * (1 – 20%). | -----   |
| <b>PCV</b>              | 5 ~ 60 cmH <sub>2</sub> O  | 1 cmH <sub>2</sub> O (based on PEEP)            | ±2 cmH <sub>2</sub> O (≤10 cmH <sub>2</sub> O); ±20% (other)  | -----   |
| <b>MV<sub>max</sub></b> | more than 18 l             | -----   | -----   | -----   |

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**

## Monitoring performance

| Item  | Range  | Resolution                    | Accuracy                              |
|---|--|-------------------------------|---------------------------------------|
| <b>V<sub>TI</sub></b>                                       | 0 ~ 2000 mL  | 1 mL                          | ±25 mL (≤100 mL);<br>±15% (other)     |
| <b>V<sub>TE</sub></b>                                       | 0 ~ 2000 mL  | 1 mL                          | ±25 mL (≤100 mL);<br>±15% (other)     |
| <b>f<sub>spont</sub></b>                                    | 0 ~ 100 bpm  | 1 bpm                         | ±2 bpm (≤20 bpm);<br>±10% (other)     |
| <b>f<sub>total</sub></b>                                    | 0 ~ 100 bpm  | 1 bpm                         | ±2 bpm (≤20 bpm);<br>±10% (other)     |
| <b>MV</b>   | 0 ~ 99 l (display value range); 0 ~ 30 l (useable) | 0.01 l (<10 l); 0.1 l (other) | ±1 l (≤3 l); ±15% (other)             |
| <b>P<sub>mean</sub>, P<sub>plat</sub>, P<sub>peak</sub></b> | 0 ~ 80 cmH <sub>2</sub> O                          | 0.1 cmH <sub>2</sub> O        | ± (2% full scale + 4% actual reading) |
| <b>PEEP</b>   | 0 ~ 30 cmH <sub>2</sub> O                          | 0.1 cmH <sub>2</sub> O        | ± (2% full scale + 4% actual reading) |
| <b>FiO<sub>2</sub></b>                                      | 15% ~ 100%   | 0.1 vol.%                     | ±3% full scale                        |
| <b>RSBI</b>   | 0 ~ 9999 1/(min × L)                               |                               | ±15%                                  |
| <b>I:E</b>  | 4:1 ~ 1:9  |                               | ±20%                                  |
| <b>Pressure waveform monitor</b>                            | -20 ~ +100 cmH <sub>2</sub> O                      | -----                         | -----                                 |
| <b>Flow waveform monitor</b>                                | -180 ~ +180 l                                      | -----                         | -----                                 |
| <b>Volume waveform monitor</b>                              | 0 ~ 1800 mL  | -----                         | -----                                 |

**CAUTION:** The monitored and displayed volume and ventilation is considered to meet precision requirement.  
The ventilator internal components which are sensitive to airflow have check valves. The replacement of the above components requires technician authorized by DIXION.

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**

## Alarm performance

| Item                                 | Range   | Accuracy   | Default value         | Remark  |
|--------------------------------------|---|--|-----------------------|---|
| <b>MV-lower limit</b>                | OFF, 0 ~ 29 l (adjustable);<br>0 ~ 20 l (useable)                 | $\pm 1$ l ( $\leq 3$ l);<br>$\pm 15\%$ (other)                                   | OFF                   | Medium priority alarm with light and buzz. Test this alarm by setting MV-lower limit different values.  |
| <b>MV-upper limit</b>                | 0 ~ 30 l (adjustable);<br>0 ~ 20 l (useable)                      | $\pm 1$ l ( $\leq 3$ l);<br>$\pm 15\%$ (other)                                   | 20 l                  | Medium priority alarm with light and buzz. Test this alarm by setting MV-upper limit different values.  |
| <b>Paw-lower limit</b>               | 0 ~ 40 cmH <sub>2</sub> O   | $\pm 1$ cmH <sub>2</sub> O ( $\leq 5$ cmH <sub>2</sub> O);<br>$\pm 20\%$ (other) | 0                     | When Paw is lower than low limit, delaying 4~15 seconds, it will occur high priority with light and buzz. Test this alarm by setting Paw-lower limit different values.  |
| <b>Paw-upper limit</b>               | 20 ~ 80 cmH <sub>2</sub> O  | $\pm 10\%$   | 40 cmH <sub>2</sub> O | When Paw is greater than high limit, it will occur high priority with light and buzz and change to expiratory state. Test this alarm by setting Paw-upper limit different values.   |
| <b>V<sub>TE</sub>-lower limit</b>    | OFF, 100 ~ 2000 mL  | $\pm 30$ ml or actual reading<br>$\pm 15\%$ , whichever is greater               | OFF                   | Medium priority alarm with light and buzz. VT low limit should be lower than VT setting value. Test this alarm by setting VTE-lower limit different values.   |
| <b>f<sub>spont</sub>-upper limit</b> | OFF, 1 ~ 60 bpm   | $\pm 2$ bpm ( $\leq 20$ bpm);<br>10% (other)                                     | OFF                   | Medium priority alarm with light and buzz. Test this alarm by setting fspont different limits.  |
| <b>T<sub>apnea</sub></b>             | 10 ~ 60 s   | $\pm 20\%$   | 20 s                  | Message prompts on display when Tapnea occurs with light and buzz. Ventilation mode switch to backup ventilation and it will switch to normal mode when next spontaneous trigger occurs. High priority alarm with light and buzz. It is validate in SIMV, PSV and SPONT. Switch to backup ventilation after alarm occurs. Test this alarm by setting Tapnea different limits. |
| <b>FiO<sub>2</sub>-lower limit</b>   | OFF, 20% ~ 99%  | $\pm 10\%$   | 20%                   | The monitored values getting by smoothing algorithm is lower than the limit for 30 seconds. Medium priority with light and buzz. Test this alarm by setting different limits.   |
| <b>FiO<sub>2</sub>-upper limit</b>   | 21% ~ 100%, OFF   | $\pm 10\%$   | 50%                   | The monitored value getting by smoothing algorithm is higher than the limit for 30 seconds. Medium priority with light and buzz. Test this alarm by setting different limits.   |
| <b>Battery Low</b>                   | Battery capability detected less than 23.4 ~ 23.7 V after startup | -----  | -----                 | Medium priority alarm with light and buzz. Test this alarm by simulation battery or battery discharge.  |

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**

| Item                                 | Range   | Accuracy | Default value | Remark   |
|--------------------------------------|---|----------|---------------|--|
| <b>Battery Exhausted</b>             | Battery capability detected less than 23 V after startup      | -----    | -----         | High priority alarm with sound and buzz. Test this alarm by simulation battery or battery discharge.   |
| <b>O<sub>2</sub> supply down</b>     | less than 0.25 MPa  | ±20%     | -----         | High priority alarm with light and buzz (in ventilation state and FiO <sub>2</sub> setting is more than 21%). Low priority alarm with sound and buzz (in standby state and FiO <sub>2</sub> setting is less than 21%). Test this alarm by adjusting oxygen source pressure.  |
| <b>Airway pressure continue high</b> | Paw is 15 cmH <sub>2</sub> O higher than PEEP for 15 seconds. | ±10%     | -----         | High priority alarm with light and buzz. Test this alarm by block exhaust vent.  |
| <b>Mains Failure</b>                 | -----   | -----    | -----         | If the ventilator without internal batteries, alarm sound keeping 120 seconds; otherwise the ventilator working by internal batteries and showing "Mains Failure" medium priority with light and buzz. Test this alarm by cutting AC power supply.   |
| <b>Leakage alarm</b>                 | -----   | ±10%     | -----         | When satisfy one of the following three conditions and delay 3 breath cycles, it will send high priority alarm with light and buzz and show circuit disconnect or leakage alarm messages.<br>1.The monitored value of inspiratory VT is 500mL higher than expiratory vale;<br>2.Paw is less than 3 cmH <sub>2</sub> O, meanwhile the monitored value of inspiratory VT is 100 mL higher than expiratory vale;<br>3.Paw is less than 3 cmH <sub>2</sub> O; meanwhile the monitored value of inspiratory VT is more than 40 mL and expiratory vale less than 5 mL. |
| <b>Alarm silence</b>                 | -----   | -----    | -----         | No more than 120 s.  |

**CAUTION:** All low limits of parameters in above table may not be set up the high limits, nor may the high limits be set below the low limits.

# PRODUCT SPECIFICATIONS VENTILATOR **Aeros-4600**

## Accessories

| Item                                | Specifications  |
|-------------------------------------|---|
| <b>Ventilation mode</b>             | Standard configuration: VCV, SIMV, SIGH, SPONT, PSV, PCV                |
| <b>Breath frequency</b>             | Standard configuration: CMV frequency 1 - 70;<br>SIMV frequency 1 - 40; |
| <b>TV</b>                           | 50 ~ 1500 ml  |
| <b>PEEP/CPAP</b>                    | Standard configuration  |
| <b>PSV level</b>                    | Standard configuration  |
| <b>FiO<sub>2</sub> adjust</b>       | Standard configuration  |
| <b>Monitor F<sub>spont</sub></b>    | Standard configuration  |
| <b>Monitor FiO<sub>2</sub></b>      | Standard configuration  |
| <b>RSBI</b>                         | Standard configuration  |
| <b>F<sub>spont</sub> high limit</b> | Standard configuration  |
| <b>FiO<sub>2</sub> high</b>         | Standard configuration  |
| <b>FiO<sub>2</sub> low</b>          | Standard configuration  |
| <b>O<sub>2</sub> Deficiency</b>     | Standard configuration  |
| <b>Battery low</b>                  | Standard configuration  |
| <b>Battery exhausted</b>            | Standard configuration  |
| <b>T<sub>p</sub></b>                | Standard configuration  |
| <b>Monitor I:E</b>                  | Standard configuration  |
| <b>Monitor PEEP</b>                 | Standard configuration  |
| <b>Alarm log</b>                    | Standard configuration  |
| <b>Screen lock</b>                  | Standard configuration  |
| <b>Alarm reset</b>                  | Standard configuration  |
| <b>Day/Night mode</b>               | Standard configuration  |
| <b>Standby</b>                      | Standard configuration  |
| <b>Unit selection</b>               | Standard configuration  |
| <b>Intended use</b>                 | Adult and pediatric   |